

Bruker AC 200 NMR Spectrometer
Command Summary for ¹H Acquisition with DISNMR Software

wallet, watch off!	C and D knobs
LINE/RGB	M
signup binder	G
"Experiment in Progress" sign face up	value [return]
spin off	M [return]
lock off	F
lift	10.5 [return]
remove reference sample	- 0.5 [return]
load positioned sample	[return]
lift off	TI [return] title [return]
spin	AZF [return] [return]
1	PXB [return] INT1.001 [return]
RJ PROTON.001 [return]	1
II [return]	[ctrl] H
O1 [return] value [return]	spin off
2	lock off
PJ PROTON.001 [return]	lift
1	remove sample
field, lock gain	load positioned reference sample
lock	lift off
lock gain	spin
Z and Z ² shims	field, lock gain
lock gain	lock
std by	lock gain
RGA [return]	Z and Z ² shims
RG [return] [return]	lock gain
ZG [return]	std by
TR [return] 2 [return]	NP [return]
EF [return]	LINE/RGB
EP [return]	"Experiment in Progress" sign face down
A and B knobs	notify next user
- / + buttons	wallet, watch, sample, spectrum
P	