

Project: Starlight Express USB Port Programming

In this project, you will need to use LabVIEW to control the tip-tilt lens from the Starlight Express. The tip-tilt has 4 motors to tilt the image in the East, West, North and South directions, and can be controlled via a serial USB port. The control is done through the USB port by sending a 7-digital string text. The following list the 4 commands to send 1-step tilt in each direction. For example, to send one step in the north:-, the text command is GN00001.

The maximum step is 45. To send 40 steps in south, the command is GS00045. You can use the LabVIEW function “**Build Text**” to build this a 7-digital string command.



The serial mode command list

Char 1	Char 2	Char 3	Char 4	Char 5	Char 6	Char 7	Result	Return
K							Find Centre	K
G	N	0	0	0	0	1	AO 1 step North	G or L (hit limit)
G	S	0	0	0	0	1	AO 1 step South	G or L (hit limit)
G	T	0	0	0	0	1	AO 1 step East	G or L (hit limit)

G	W	0	0	0	0	1	AO 1 step West	G or L (hit limit)
---	---	---	---	---	---	---	----------------	-----------------------

Write a code to control each direction. In this assignment, you need to use the random function to generate a random value to control the amplitude in each direction that can be chosen and changed from the front panel. A Led light will switch on to show the current direction being chosen. This can be done by using a while loop. See the following picture for the LabVIEW Front Panel.

