

CSUN BS, Radiologic Sciences Program
Clinical Competency Assessment for General Specials
 (Use one form per competency)

The student must pass all sections of this form. In every section, each student must accumulate the minimum points necessary based on their Junior or Senior/Post Grad status. If a licensed RT intervenes during any part of the evaluation, the assessment is discontinued and the student must be re-evaluated on a different exam. The following total points are needed to pass this competency: Juniors **52** and Seniors/Post Grads **59**.

Name: _____ Procedure: _____

Date of Exam: _____ Tech's Name and Exam Time: _____

Scale: 0=Unacceptable; 1=Needs Improvement; 2=Competent/Acceptable; 3=Above Average; 4=Exceeds Expectations (at RT level)

I. Patient Care & Communication					
1. Verifies correct patient using two patient identifiers, procedure, and clinical diagnosis/indication.	0	1	2	3	4
2. Communication skills – Provides clear and complete explanation of procedure to age appropriate patient, surname use, etc.	0	1	2	3	4
3. Obtains patient assessment, preparation, history, and consent. (NPO, labs, contrast allergies, etc)	0	1	2	3	4
4. Insures patient privacy and dignity.	0	1	2	3	4
5. Knows emergency protocol for Code situations.	0	1	2	3	4
6. Insures patient safety and comfort throughout procedure.	0	1	2	3	4
7. Practices safe patient care methods including falling precautions.	0	1	2	3	4
8. Knows how to operate bedside controls and remembers to put side rail up after the exam.	0	1	2	3	4
9. Provides proper breathing instructions; appropriate for exam and patient language.	0	1	2	3	4
10. Preps area of interest using proper sterile technique and maintains the sterility of the field throughout the procedure.	0	1	2	3	4
Total for Section I: _____					
Maximum points = 40 Minimum passing points: Juniors=16 and Seniors/Post Grads=18					
II. Technical Requirements & Positioning					
1. Correctly enters patient information.	0	1	2	3	4
2. Selects appropriate exam protocol or adjusts as necessary.	0	1	2	3	4
3. Adjusts exposure factors and selects proper image receptor.	0	1	2	3	4
4. Positions patient efficiently and correctly for the examination. (Scout, overheads, etc.)	0	1	2	3	4
5. Uses proper immobilization devices.	0	1	2	3	4
6. Uses tube angulations properly.	0	1	2	3	4
7. Knows storage location of necessary supplies.	0	1	2	3	4
8. Examination completed in a timely fashion.	0	1	2	3	4
9. Properly archives images to PACS or prints images.	0	1	2	3	4
10. Uses proper SID.	0	1	2	3	4
11. Uses correct image sizes.	0	1	2	3	4
12. Uses proper alignment of x-ray tube to image receptor.	0	1	2	3	4
13. Uses proper central ray location.	0	1	2	3	4
14. Uses proper orientation of Image Receptor. (Lengthwise, crosswise)	0	1	2	3	4
15. Properly selects and prepares sterile procedure tray.	0	1	2	3	4
Total for Section II: _____					
Maximum points = 60 Minimum passing points: Juniors=24 and Seniors/Post Grads=27					
III. Image Critique					
1. Diagnostic image quality (contrast/density), exposure range must meet department protocols.	0	1	2	3	4
2. Patient positioning demonstrates anatomy and pathology properly.	0	1	2	3	4
3. Able to identify which projection of the image was taken.	0	1	2	3	4
Total for Section III: _____					
Maximum points = 12 Minimum passing points: Juniors=5 and Seniors/Post Grads=6					
IV. Radiation Protection					
1. Patient assessment. (Pregnancy, child bearing age, adolescent, etc.)	0	1	2	3	4
2. Uses proper shielding and radiation protection for patient\staff\parents.	0	1	2	3	4
3. Uses proper collimation of the body part.	0	1	2	3	4
4. Documents fluoroscopy time and fluoroscopy air kerma/DAP appropriately.	0	1	2	3	4
Total for Section IV: _____					
Maximum points = 16 Minimum passing points: Juniors=7 and Seniors/Post Grads=8					

***** If the student receives a "0" in any of the above categories, then he/she will not receive credit for the comp. *****

RT's Name (Print): _____ RT's Signature: _____ Date: _____

In order for this form to be valid, the technologist must be registered by the ARRT for at least 2 years.